

**SCHILLER**

www.schiller.ch

CARDIOLOGY**BIO SERVE PROMOTION NEW**

FRED PA-1 (Public Access 1)

Dedicated to public access, SCHILLER's latest defibrillator, the FRED PA-1, has been designed so that even untrained users can save lives. By lifting up the device cover, the FRED PA-1 starts up immediately and guides the rescuer step-by-step during the entire resuscitation process.

**Code: 1-127-9902****FRED PA-1****Automatic/ Semi-Automatic Defibrillator.**

Public Access Defibrillator with Biphasic Truncated Exponential waveform (BTE), metronome and validity Radio-Frequency Identification (RFID) detector.

Includes following accessories:

- LiMnO₂ Lithium non-rechargeable battery
- 1 set of pre-connected adult pads (disposable)

**10 year warranty
& 6 years
battery life!**

Features:

- 10 year warranty period & 6 years battery life with weekly self checks
- Clear user interface: 1-2-3 steps with spoken instructions and clear pictures
- Pre-connected electrodes for faster application to the patient's chest
- Automatic self-tests for detection of electrode expiration and battery capacity with innovative Radio-Frequency Identification (RFID) Technology
- The FRED PA-1 is available either as;
 - ⇒ Automatic (delivers the shock without any action of the rescuer)
 - ⇒ Semi-Automatic (prompts the rescuer to deliver the shock by pressing the shock button)
- Easy to carry hard case
- Data Management: up to 2 hours of intervention including 500 events can be stored in a simple SD card.
- For all types of patients with adult & children electrodes
- BTE: Biphasic Truncated Exponential waveform
- AED energy: 150/ 200/ 200J for adults
50/ 50/ 50J for children
- Configurable metronome pace for CPR

PRICE ► \$2, 995*

(09) 846 4802 **(09) 846 6808** **sales@bioserve.co.nz** **www.bioserve.co.nz**

BIO SERVE
Medical Sales

100% NZ
Owned &
Operated!

Bio Serve NZ Ltd
PO Box 56-516, Dominion Rd, Auckland 1446,
31 Morningside Dr, Mt Albert, Auckland 1025, NZ.

***Prices Freight and GST Excl.**

Pub. 12-Oct-16 **Valid until 31/12/2016**