



# *Waterlase*® Tips Quick Reference Guide



The following guide lists the most widely recommended BIOLASE single-use tips by procedure type. The purpose of this guide is to inform the clinician of our most highly utilized single-use tips by procedure category to help guide which tips may be better for certain patients. All patients are different and recommendations may vary, so this is meant to be utilized as a guide versus protocol.

Contact your **BIOLASE Distributor** for more information or to check out the latest offerings.



*Note: Images enlarged to show detail.*

## RESTORATIVE

### CLASS I-II

MZ6 – 6mm

MZ8 – 6mm

### CLASS III - V

MZ5 – 6mm

MZ6 – 3mm

MZ6 – 6mm

### DECIDUOUS

MZ5 – 3mm

MZ6 – 3mm



*MZ5– 6mm*

*Longer tips may be required depending on access. Contact your BIOLASE distributor for additional sizes.*

## SOFT-TISSUE

Frenectomy, biopsy, gingivectomy, gingivoplasty, troughing, and hemostasis

MZ5 – 6mm

MZ6 – 6mm



*MZ6 – 6mm*

*Longer tips may be required depending on access. Contact your BIOLASE distributor for additional sizes.*

## PERIO TIPS

REPAIR™ Perio Protocol, osseous crown lengthening (open or closed), and flap surgery for removal of pathology

MZ6 – 6mm


MZ6 – 14mm

RFPT5 – 14mm



*RFPT5 – 14mm*

*Longer tips may be required depending on access. Contact your BIOLASE distributor for additional sizes.*



*Radial Firing  
Perio Tip*

## IMPLANT

REPAIR™ Implant Protocol, implant uncovering, socket debridement, and flap surgery for removal of pathology

MZ6 – 6mm

MZ6 – 14mm


RFPT5 – 14mm

SFT8 – 18mm (“Side Firing Tip”)\*

*\*Side Firing Tip is a premium, sapphire tip that is reusable versus single-use.*

*Longer tips may be required depending on access.*

*Contact your BIOLASE distributor for additional sizes.*





## ENDO

Pulp chamber access, pulpotomy, root canal cleaning and shaping

MZ6 – 9mm

RFT2 – 25mm



Contact your BIOLASE distributor  
to order tips today!



4 Cromwell, Irvine, CA 92618 USA  
biolase.com • +1 949.361.1200  
international@biolase.com • NASDAQ: BIOL

MADE IN THE USA

©2019 BIOLASE, Inc. All rights reserved. 17-0946